

# General Purpose Multi-Tissue Ultrasound Phantom

*The standard for ultrasound quality assurance  
Two Phantoms in one case*

The CIRS series of ultrasound phantoms, unlike human subjects or random scannable materials, offers a reliable medium which contains specific, known test objects for repeatable qualitative assessment of ultrasound scanner performance over time.

This phantom is constructed from the patented solid elastic material, Zerdine™.<sup>(1)</sup> Zerdine™, unlike other phantom materials on the market, is not affected by changes in temperature. It can be subjected to boiling or freezing conditions without sustaining significant damage. Zerdine™ is also more elastic than other materials and allows more pressure to be applied to the scanning surface without subsequent damage to the material. At normal or room temperatures the Zerdine™ material found in the Model 40 will accurately simulate the ultrasound



Model 040

---

## Complies with AIUM Standard for Quality Assurance.

---

characteristics found in human liver tissue. The Model 40 was designed to allow for assessment of uniformity, axial and lateral

resolution, depth calibration, dead zone measurement, and registration within two different backgrounds of 0.5 and 0.7 dB/cm/MHz.

<sup>©</sup>US PAT# 5196343

**CIRS**  
Tissue Simulation Technology

2428 Alameda Avenue • Suite 212 • Norfolk, Virginia 23513 • USA

# Model 40 Specifications:

**MATERIAL:** Zerdine™(1), solid elastic water-based polymer  
 Freezing Point: 0° C  
 Melting Point: Above 100° C

**ATTENUATION COEFFICIENT:**  
 0.5 dB/cm/MHz  
 0.7 dB/cm/MHz

**SPEED OF SOUND:**  
 1540 m/s

**SCANNING WELL:**  
 1cm deep

**SCANNING MEMBRANE:**  
 Saran

**TARGETS:**  
 Material: Nylon  
 Monofilament  
 Wire Diameter: 0.1mm

**VERTICAL PLANE TARGETS**  
 Number of Groups: 1  
 Number of Targets Per Group: 16  
 Depth Range: 18cm  
 Spacing: 1cm


**HORIZONTAL PLANE TARGETS**  
 Number of Groups: 2  
 Number of Targets: 4 and 7  
 Depth Range: 3cm and 9cm  
 Spacing: 1cm and 2cm

**RESOLUTION TARGETS:**  
 Number of Arrays: 4  
 Depths: 2.5cm, 6cm and 10cm  
 Axial Intervals: 0.5, 1, 2, 3, 4, and 5mm  
 Horizontal Intervals: 1, 2, 3, 4, and 5mm

**CYSTIC TARGETS:**  
 Number of Targets: 4  
 Diameter of Targets: 2, 4, 6 and 8mm  
 Depth of Targets: 2, 4, 6, and 8cm  
 Attenuation: <0.07dB/cm/MHz  
 Speed: 1540 m/s  
 Contrast: anechoic

**HIGH CONTRAST TARGETS:**  
 Number of Targets: 4  
 Diameter of Targets: 2, 4, 6 and 8mm  
 Depth of Targets 2, 4, 6, and 8cm  
 Attenuation: 1.0 dB/cm/MHz  
 Speed: 1540 m/s  
 Contrast: +15 dB v.s. background

**HIGH DENSITY TARGET:**  
 Material: PMMA  
 Diameter: 1/16"  
 Depth: 6cm



Phantom comes with detachable scanning wells to accommodate large sector probes and small endocavity probes. It is packaged in a hermetically sealed, air tight, rugged carrying case.